
Aseptik-Klemmverbindungen DIN 11864-3, Form A

Aseptic clamp pipe connections DIN 11864-3, Form A

Hygiene-Klemmverbindungen DIN 11853-3

Hygienic clamp pipe connections DIN 11853-3

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Application

The aseptic and the hygienic clamp pipe connections are compatible and differ only in the axial length. They are used with process facilities for liquid media with high purity requirements. Preferential fields of application are bio-engineering, cosmetics, pharmaceutical, chemical, food, and beverage industry.

Assembly

During the assembly the O-Ring is positioned in the clamp ferrule with groove whereby the O-Ring is automatically held. Then the clamp ferrule with notch is put against so that both ferrules interlock. By this the ferrules and the O-Ring are coaxially centred. Then the clamping ring is applied to the flanges and closed (alternatively with safety catch) until the fronts of the flanges get together. Thereby, the O-Ring is deformed in such way that an even and free of dead space thoroughfare and a tight connection is attained.

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Articles	Page
Clamping rings	196
O-Rings	196
Clamp blind flanges with groove	200
Clamp blind flanges with notch	201
Sight glasses	202
Clamp ferrules with groove (NS)	203
Clamp ferrules with notch (BS)	205
Intermediate pieces	207